

SOP TITLE: CONFIGURATION OF OFFLINE CHT APPLICATION SETTINGS IN KENYAEMR

Objectives

This document should help you to quickly configure CHT Application settings inside KenyaEMR

Target audience

Facility champion, KenyaEMR Technical lead, Palladium and other stakeholders

Applicability

This procedural document is to be used by Facilities that integrating KenyaEMR and CHT Medic Mobile Application

Prerequisites:

- CHT App must be installed
- A user with administrative role must be setup in the CHT App under the Users inside App Management

Step 1:

Log into KenyaEMR with the administrator account:

Step 2:

Click the home tab:

Step 3:

Click the developer tab:



Step 4:

Click the **Legacy admin UI** button:



Step 5:

Go the **Maintenance** section:

Maintenance

Set Implementation Id
System Information
View Quick Reports
Settings
Advanced Settings
View Server Log
View Database Changes
Manage Locales And Themes
View Logged In Users
Search Index

Step 6:

Click the **Advanced Settings** section:

Advanced Settings

Step 7:

Locate the following settings in the application

/h	
cb6f4d4b-73cc-4c42-97cb-0db5a631190a	Ren
https://cht.hmislocal.org/medic/_bulk_docs	Ren
	li.
lk.	
medic	Ren
4	
	https://cht.hmislocal.org/medic/_bulk_docs

Kindly set the following settings based on the user with administrative rights you created for the CHT application installed.


Medic.chtUser:-This implies the username based on the user with the administrative rights you created in CHT app.

<u>Medic.chtPwd</u>:-This implies the password based on the user with the administrative rights you created in the CHT App.

Medic.chtServerUrl: This should point to the IP address where the CHT application is running. Kindly change the localhost to the IP address where CHT app is running if it is installed in a different machine.

CONFIGURATION NEED IN THE CHT APPLICATION


Step 1: Login in the medic mobile with the administrative password in this case it is: UserName: medic Password: cb6f4d4b-73cc-4c42-97cb-0db5a631190a



Step2

Access the fauxton user interface which is a native web-based interface built into CouchDB with changing the url to:

https://cht.hmislocal.org/ utils




Step 3

Go and click the configuration tab:



Step 4

Click the Add Option tab and add the following options as shown below:



The Options needed to be added is:

Section: medic-credentials

Name: kenyaemr Value:Admin123

The final option is to click the Create button that adds the section ,option and value as part of the configuration.

